

# Image-Pro® Insight

Image Analysis Made Easy

Capture, Process, Measure, and Share

# Image-Pro<sup>®</sup> Insight

**Image Analysis Made Easy** 


#### Capture, Process, Measure, and Share

Image-Pro Insight, the latest image analysis software from Media Cybernetics, offers a wide range of tools for capturing and analyzing images.

Unlike most imaging systems, which require extensive training, Image-Pro Insight was designed for quick startup and ease of use. Utilizing more than 25 years of user input, Image-Pro Insight offers intuitive tools that make it easy to capture and analyze images and extract valuable data from them.

# Capture with Ease

Capturing high-quality images is made easy with a range of acquisition tools. Image-Pro Insight supports an ever-growing list of scientific cameras along with TWAIN (32-bit) acquisition to ensure a superior out-of-the-box capture experience.



### Saturation Warning on Live Preview

Receive an automatic visual cue to prevent you from acquiring saturated images.



#### Live Extended Depth of Field – No Automation Needed


Capture beautifully focused images, regardless of depth, with Live EDF. No automated z motor required.





## Apply High Dynamic Range (HDR) during Capture

More faithfully capture the widest range of intensity present during image acquisition.




#### Live Tiling – No Automation Needed


Stitch together multiple fields of view as you manually move around your specimen with Live Tiling. Combine with Live EDF to create a focused, tiled image. No automated stage required.

# Process & Analyze

Implement sophisticated image processing and enhancement functions with two clicks or less. Image-Pro Insight offers a wide range of tools that make it easy to analyze and reveal important details within your images.



Analyze with Line Profile Add multiple line profiles to live or captured images.




Create Color Composite Images Build composite color images by simply right-clicking on each grayscale image to select a color tint.






# Auto-align and Apply Extended Depth of Field

Correct for microscope shift with post-acquisition auto-alignment tools. Extended Depth of Focus (EDF) creates a focused image from a series of images captured at different z planes.




# Extract & Merge Color Channels

Easily combine and separate images into RGB, HSI, and HSL color channels.



# View Image Intensity with Surface Plot

Display a 3D graph of image pixel intensity with the Surface Plot tool. Visualize using texture, color, and wireframe settings, and create instant movies.



#### Correct for Noise & Enhance Detail with Filters

Remove background noise and reveal hidden details with an extensive set of image processing filters. Preview filters on your active image for instant results.




#### Navigate & Edit Time-Lapse Images in Gallery View

Navigate through time-lapse and image sequence files with the Gallery View tools.


# Measure & Quantify

Extract accurate data from your images with a wide range of measurement tools.




#### **Classify Measurements**

Use different colors, shapes, and custom names to classify and organize your measurement data.




### Apply Semi-automated Measurements

Segment, split, join, shrink and grow objects with just a couple of clicks using the Magic Wand and the Relative Measurement tools.



#### **Apply Caliper Measurements**

Use the Line Profile tool to measure pixel intensity along a line.



### Export Measurement Data to Excel

Export measurements and data graphs to Excel.


# Report & Share

With Image-Pro Insight, it's easy to present your images and analysis data. Use one-click tools to export images, measurements, and annotations to Microsoft® Excel® as well as to popular presentation programs like PDF, Microsoft Word®, and Microsoft PowerPoint®.



### Create Reports in PDF, Word & Excel

Share your images and data in easy-to-use report templates.



#### **Annotate Images**

Point out features of interest with customizable annotation tools. Easily adjust fonts, colors, shapes, and more. Save these overlays on your images to share with others.

### Join the Image-Pro® **Analysis Community**

#### **Get Started with Free Video Tutorials**

As an Image-Pro Insight user, you'll have access to a growing library of free 3- to 5-minute video tutorials. Learn timesaving tips and familiarize yourself with the wide range of analysis tools available to you – without ever having to leave your computer.

#### **Technical Support**

Receive installation and configuration troubleshooting from our technical support engineers via phone or email. Purchase an Image-Pro Insight Maintenance Agreement to receive a full year of technical support and software updates.

#### Image-Pro Users Email List

When you choose Image-Pro software, you instantly join a community of image analysis professionals. Join the free Image-Pro users email list to ask questions and learn from other users.

#### **Upgrade to Image-Pro Premier**

If you find your analysis needs growing beyond the tools in Image-Pro Insight, you can easily upgrade to Image-Pro Plus, which offers macro customization, automated counting/ sorting, and other advanced options. Contact your local reseller or Media Cybernetics for details.

# Free Video **Tutorials**

#### **Designed To Meet Your Requirements**

#### Life Science Applications



**Animal Behavior** 

**Food Science** 

Marine Biology

#### **Industrial Applications**











#### **Features Overview**

#### Capture

Single frame & movie acquisition

Live tiling & EDF

Crosshair grid on live image display

TWAIN acquisition (32-bit)

**HDR** 

#### Process & Analyze

Post-acquisition alignment & EDF

Frame averaging

Color composite

Pseudocolor

Filters – enhancement, edge, large, morphological

Binary mask

Sequence / movie editing

#### Measure & Quantify

Advanced line profile analysis

Semi-automated measurements

Measure live images

Snap measurements

Best-fit measurements

Classify measurements

Interactive data table

Data histogram

#### Report & Share

Annotation tools

Quick save for publication & analysis

Export to Excel, PowerPoint & Word

Reporting tools (including PDF)

Create snapshot of zoomed view

Create snapshot with overlays

#### View

Surface plot view & movie creator

Gallery view

Lock & zoom multiple images

Zoomed image navigator

Full-screen mode

#### **System Requirements**

#### **Supported Operating Systems:**

Windows® XP SP 3 and later – 32-bit

Windows® 7 Professional and Ultimate- 32 & 64-bit

Windows® 8 Professional and Enterprise-32 & 64-bit

#### **Minimum Requirements**

- Dual-core processor or better
- 4 GB of RAM
- 5 GB available on installation drive + free space for images
- DVD-ROM drive
- 1,024x768 display with 32-bit video card
- USB port required for hard licenses and offline license activation
- Internet connection required for online services\*
- Internet Explorer version 9 or higher

#### **Recommended Requirements**

- Quad-core processor or better
- 16 GB of RAM memory or higher
- 20 GB available on installation drive + free space for images
- DVD-ROM drive
- 1,280x800 display with 1GB video card
- USB port required for hard licenses and offline license activation
- Internet connection required for online services\*
- Internet Explorer version 9 or higher

Microsoft .NET Framework 4.0 (automatically installed by application).

Specifications are subject to change.

#### **Supported File Formats**

Description	Extension(s)	Write	Read
TIFF	*.tif, *.tiff	X	X
JPEG	*.jpeg, *.jpg	X	Х
JPEG 2000	*.jp2, *.jpf	X	×
BBMP	*.bmp	X	Х
PNG	*.png	X	Х
Flex Image Transport System	*.fts, *.fts, *.fit, *.fits	X	×
Image-Pro Set	*.ips	X	Х
Image-Pro Raw	*.ipraw		×
Image-Pro Sequence Files	*.seq		×
AVI	*.avi	X	×
Windows Media Video	*.wmv	X	Х
MPEG 4	*.mp4, *.mov, *.m4v	X	X
BIGTIFF	*.btf	X	×



Media Cybernetics, Inc. North & South America: + 1 301 495 3305 Europe, Middle East & Africa: +44 (0)1628 477025 China: +(86)21 33773539

Asia, Australia, India & New Zealand: + 65 6408 6245

Japan: +81 3 5639 2751

www.mediacy.com | info@mediacy.com

Copyright © 2014 Media Cybernetics, Inc.
All rights reserved.
Media Cybernetics, the Media Cybernetics
logo, and Image-Pro are registered
trademarks of Media Cybernetics, Inc.
All other trademarks are used with
permission and remain the property of
their respective owners.

<sup>\*</sup>This product may allow you to access certain features that are hosted online ("online services"), provided you have a high-speed Internet connection. The online services include but are not limited to, automatic updates, support links, access to video tutorials, and more