

Comparison of CSU series

Model		CSU-X1		CSU22	CSU10
		High-end	Basic		
Imaging Speed (Max. fps)		2,000	360	1,000	360
Scanner Motor Rotation Speed (rpm)		1,500-10,000 (Variable)	1,800 (Fixed)*1	1,500-5,000 (Variable)	1,800 (Fixed)
Recommended camera exposure time		0.5msec <	33msec <	1msec <	33msec <
Rotation position trigger signal		External signal output possible	None*1	External signal output possible	None
Filters	EX	Option	Option	Supplied up to 3 filters	Supplied 1 filter
	DM	Option (up to 3 filters)	Option (1 filter)	Supplied up to 3 filters	Supplied 1 filter
	EM	Option (up to 6 filters with filter wheel)	Option (1 filter)	Supplied up to 3 filters	Supplied 1 filter
Addition or exchange of filters		At user site : DM block and filters (EX, EM.) At Yokogawa factory : DM		At Yokogawa factory	At user site
*1 option					